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INSTRUCTIONS FOR PROFESSIONAL USE IN THE HEALTH / MEDICAL / DENTAL FIELD FACE MASK FOR PROTECTION AGAINST COVID-19 - PPE-R/02.075-V2

This filtering half mask is manufactured for COVID-19 protection. Suitable for protection in medical, dental, sanitary applications even in presence of aerosol. This filtering half mask is not a filtering half mask for general use (chemicals, asbestos, etc.) and shall not be used for purposes other than protection against COVID-19. Warning - this product is not flame resistant and must not be used in areas with open flames. AXELMED FPS SAFECOMFORT is a medical device that meets the requirements of EN 14683 Type IIR, filtering face masks against bacteria, viruses in air dispersion and in aerosol.

Pay particular attention to warning phrases where indicated. WARNINGS and LIMITATIONS Always make sure that the product is: Suitable to be used; Worn correctly; Worn for the entire period of possible exposure to the infection; Replaced when needed; For the proper use of the device, the user should be familiar with the correct methods of use in order to protect against certain airborne contaminants. Failure to comply with the instructions on the use of the device to protect the respiratory tract and / or the failure to wear for the entire period of exposure can cause serious damage to health and disability, including permanent ones. For the correct use of the device, it is advisable to follow local regulations and refer to all the information provided. For further information, contact an expert at work safety or an Axelmed representative; Before use, the user must be trained in the correct use of the product; in accordance with the standards and guidelines in force on health and safety; In these products there are no components produced in natural latex; These products do not protect against gas and vapors. The use in atmospheres containing less than 20% oxygen is not recommended; Do not use to protect the respiratory tract against unknown or immediately life-threatening atmospheric contaminants; **ATTENTION:** beard, sideburns or mustache could prevent a good seal of the mask on the face. Leave the potentially contaminated area immediately if breathing becomes difficult or dizziness or other discomfort appears; Eliminate and replace the face mask if it is damaged or if the respiratory resistance becomes excessive; Do not modify or repair the face mask; Dispose of the mask in accordance with local regulations; Before use, always check that the product has not expired.

LEGEND SYMBOLS CEN/ISO

	Lot Number		DO NOT use after expiry date		Use only if packaging is not broken
	Article Code		NOT Sterile		Caution: read instruction before use
	Manufacturer		Do not expose to sun rays		MAX 90%

STORAGE AND TRANSPORT Store the product following the instructions provided by the manufacturer inside the package. Average conditions can exceed 25°C for limited periods. Use the original packaging to store and transport this product. **Do not store in direct sunlight.** The expiry date is in the format: MM / YY. Dispose of in accordance with local regulations.

TECHNICAL FEATURES MEDICAL DEVICE CE EN 14683 TYPE II R, REGISTERED IN THE MINISTRY OF HEALTH (IT), BACTERIAL FILTRATION EFFICIENCY >99.9%, WITHOUT EXPIRATION VALVE, AUTOCALVABILE, DEVICE SUITABLE FOR USE IN COVID-19 INTENSIVE CARE UNIT, SURGERY, MEDICAL AND DENTAL PRACTICE, EMERGENCY AND AMBULANCE, RSA; Professional surgical protection mask entirely manufactured in Italy by AXELMED; 5-layer laminated material in pure TNT PP made in Germany for AXELMED. The hypoallergenic does not absorb halitosis and exhaled humidity. Conformity with EN 14683 TYPE IIR (R means resistant to blood splashes antislava); Effective barrier against viruses, bacteria and spores in air dispersion and aerosols; BFE 99.9% bacterial filtration efficiency tested with *Staphylococcus Aureus* with an average diameter of 0.7 - 1.2 microns; Excellent breathability with spongy intermediate layer to dissipate exhaled heat; Reusable for months and machine washable (40°C - 35 minutes); ATTENTION: due to the anti-splashes property of the outer layer, the stamp on the mask may fade after washing. This does not alter the filtering properties; One minute disinfection with a domestic microwave; The innovative SAFECOMFORT design adheres perfectly and does not cause irritations Soft elastic tubulars adjustable with the self-adaptive sliding system (very long; can be tightened outside the hood of the protective suit, protective visors or headgear) Absence of exhalation valve to AVOID virus transmission even when worn by positive people COVID-19 or patients at risk. The special shape of the nose piece allows you to place the protective glasses on the mask to prevent fogging. Bacterial filtration efficiency (IN and OUT) BFE 99.9% >99.9% after 5 cycles in the washing machine 40°C x 60' - >99.9% after 20 cycles in autoclave Class B 121°C x 28'. Breathability <60 Pa / cm (53.9 Pa / cm² after 5 cycles in the washing machine 40°C x 60' - 55.1 Pa / cm² after 20 cycles in class B autoclave 121°C x 28'). Bioburden: <30 cfu / g, exceed ISO EN 11737 (bacterial contamination after production; 12 cfu / g - Yeast 2.6 cfu / g - Yeast <0.4 cfu / g). Biocompatibility meets the requirements of the ISO 10993-1 face mask. Splash Test >16.0 kPa according to EN 14683 Type IIR.

USE AND MAINTENANCE INSTRUCTIONS In order for the AXELMED FPS SAFECOMFORT protection mask to be effective it is necessary to use it in combination with frequent hand cleaning performed with the help of disinfectant detergents, it is necessary that it is not used by several people and it is necessary to use it correctly. Specifically, it is advisable to follow the following steps: clean your hands with an alcohol or soap-based detergent and water before putting it on; when using for the first time, it is possible to adjust the two right and left cotton pads in the nose pad by pressing them or pulling them out according to the shape of the face and, if necessary, in special cases (small people, children), also tighten the width of the elastic with a knot; under the lower part of the elastic and insert it around the neck, then position the upper part of the elastic on the nose; the elastic is sliding and adapts independently even if you wear a face protection visor and the headress of the biological protection suit; cover the mouth and nose with the mask making sure that there are no spaces between the face and the mask along the entire perimeter of contact; it is recommended to keep the mask as high as possible and close to the eyes for perfect closure; avoid touching the mask externally while using it (if you do, clean your hands with an alcohol or soap and water based detergent); to remove the mask, free the upper elastic from the nose towards the front and then widen the sliding part behind the neck without touching the external part of the mask.

VERIFICATION OF THE SEAL cover the front of the face mask with both hands making sure not to compromise the seal on the face; breathe out quickly; if air leaks are felt around the nose, reposition the face mask higher to eliminate the leak. Repeat the leak test as above; if air exits the edges of the respirator, adjust the elastic bands on the sides and on the head to eliminate the leak. Repeat the leak test as above; in the event that a perfect seal cannot be obtained, DO NOT enter into the contaminated area. The face seal of the face mask should be checked according to the national legislation in force, if any.

MAINTENANCE OF AXELMED FPS SAFECOMFORT WASH it is recommended to keep the mask as clean as possible. If it gets dirty externally after several uses it can be washed in the washing machine by inserting it in a tightly closed white cotton sachet or liner, for a maximum of 5 washes, with the aid of mild washing machine detergents without bleach, possibly also disinfectants (rapid cycle max 35° max 40°C). Before washing, remove the two cotton swabs from the hake. Hanging drying is recommended, at room temperature and away from heat sources and sunlight.

ATTENTION: Washing with different settings can cause deformations and decrease the filtering effectiveness. ATTENTION: due to the water-repellent properties of the outer layer, the stamp on the face mask may fade after washing. This does not affect the indicated filtration values.

DISINFESTAZIONE / STERILIZZAZIONE A major limitation of protective face masks against external pathogens is the fact that, for hygienic reasons, masks are mostly disposable medical devices. This entails a series of problems with which to face such as, for example, the difficult availability of the product in conditions of strong demand from the market, an increase in costs, greater environmental pollution. These reasons have led to discussions recently about possible procedures useful for disinfecting / sterilizing masks at high temperatures and the Institute of Health is working hard to agree on a safe process to be able to reuse them. To date, however, there is no standard and accepted procedure for disinfecting surgical masks and facial filters (FFP3) and, very often, the quality or type of materials used for the realization of the devices do not allow to use disinfection / sterilization methods without altering their shape and / or function. In order to evaluate the behavior of the AXELMED SAFE COMFORT facial mask in conjunction with a sterilization process, it was decided to test the product by using a class B autoclave. The device underwent 20 sterilization cycles at 121°C for 28 minutes and tests were then performed to evaluate BFE with results greater than 99.9%. By summarizing the disinfection / sterilization methods to be adopted, the following are to be taken into consideration according to the various areas: Surgical / hospital / medical field: sterilization in heat-sealed pouch with standard fractionated vacuum cycle (MAX 121°C - MIN 8' / MAX 28') in Class B autoclave; First aid / ambulance / medical guard: it can be sprayed with preferably volatile non-detergent alcohol-based disinfectant spray (eg isopropanol) and left to dry completely without touching it, open and at room temperature.

REUSE The AXELMED FPS SAFECOMFORT mask, although falling under the ISO 14683: 2019 standard, unlike the disposable surgical masks for which a single use is recommended for each patient treated; however not exceeding one work shift without removing them with a maximum of 4 hours, it is instead a medical device CE REUSABLE therefore wearable for longer periods: all day and up to over a month of daily use. In order to reuse AXELMED FPS SAFECOMFORT surgical masks, it is necessary to carry out disinfection / sterilization procedures. Unlike what happens with many surgical masks or PPE which, once subjected to sterilization or disinfection cycles, lose part of their filtering capacity, in case of AXELMED FPS SAFECOMFORT the tests performed confirm that even after 20 sterilization cycles at 121°C the filtering capacity remains unchanged and equal to 99.9%. From the point of view of filtering capacity, the tests carried out (before and after the sterilization cycles) found a BACTERIAL FILTRATION of 99.9%. Wanting to summarize what has just been explained, in the light of the tests carried out and the characteristics of AXELMED FPS SAFECOMFORT, it is possible to identify the device as a surgical mask (being EN ISO 14683 certified) with the possibility of using more than one PPE (being reusable several times and for longer periods of time).

APPROVALS These products meet the safety requirements of the European Directive and therefore carry the CE 2017/745 marking. Version and emission data of these Operating Instructions: V03-EN-09/2020 - AXELMED S.R.L. is a Company Certified ISO 13485: 2016 FFS Model registered in the EU (EUIPO) No. 007776133-001 Patent pending (PCT). Declaration of Conformity online at www.axelmed.com

AXELMED S.R.L. is certified ISO 13485: 2016 with Cert. N° S/0020-2016 and registered as a manufacturer of Class I, IIa, IIb Sterile Medical Devices at the IT Ministry of Health. Quality Assurance System compliant with the requirements of Annex II of Directive 93/42/CEE amended by Directive 2007/47/CE and implemented with Legislative Decree 46/97 and subsequent amendments.

Scheda Tecnica - Technical Data Sheet

SAFE COMFORT

FP5 SAFECOMFORT

PPE FACE MASK FOR PROTECTION AGAINST COVID-19 (PPE-R/02.075-V2)

Description

- Material: Non Woven Fabric, Polypropylene, 5 Layers with different density/thickness/texture
- Color: White - Packaging: 1 p/box - MADE IN ITALY
- CE Medical device registered at the Italian Ministry of Health with No. 1944143 / 1948544 / R
- EU Registered Model EUIPO No. 007776133-001 - (PCT PL.P.)

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Version 003/SP2020 © 2020 AXELMED SRL - All rights reserved. Reproductions, translations, storage and electronic processing are not permitted without prior written permission from AXELMED SRL. Changes, errors and misprints under reserve. As products and regulations are under continuous development all illustrations, texts, data, features, etc. may vary without notice.

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AXELMED S.r.l.