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**Novel Coronavirus Antigen Detection Kit (Colloidal Gold)**  
**Ref.: COVID-19-NG02**  
**Version: 04**  
**Specifications: Sputum/Stool**  
**Effective Date: 2020-10**

**INSTRUCTION FOR USE**

**For professional use in vitro diagnostic use only**

**PRODUCT NAME**  
 Novel Coronavirus Antigen Detection Kit (Colloidal Gold)

**PACKING**  
 1 package, 25 pieces/box

**INTENDED USE**  
 This product is suitable for the qualitative detection of novel coronavirus in stool or sputum sample. It provides an aid in diagnosis of infection with novel coronavirus.



**SUMMARY**  
 The novel coronavirus belonging to the beta genus, COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection. Asymptomatic, infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough, nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

**PRINCIPLE**  
 The Novel Coronavirus Antigen Detection Kit (Colloidal Gold) is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid protein from SARS-CoV-2 in throat swabs. The test strip is composed of the following parts: memory sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal gold conjugated with the monoclonal antibodies against the nucleocapsid protein of SARS-CoV-2. The reaction membrane contains the secondary antibodies form nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen presents in the sample, a complex formed between the anti-SARS-2 conjugate and the virus will be captured by the specific anti-SARS-2 monoclonal antibodies coated on the test line region (T). Absence of the C line suggests negative result. To serve as a procedural control a read line will always appear in the control line region (C), indicating that proper volume of sample has been added and membrane wicking has occurred.

**COMPOSITION**  
 Each pack of rapid test kit includes:  
 1) Disposable test card  
 1) Disposable sample extraction tube  
 1) Cotton swab  
 1) Disposable paper cup


**STORAGE AND STABILITY**  
 1. Store as packaged in the hermetic bag at the temperature (2-30°C or 38-65°F) and avoid direct sunshine. The kit is stable within the expiration date printed on the labeling.  
 2. Once open, the hermetic bag, the test should be used within one hour. Prolonged exposure to hot and humid environment will cause product deterioration.  
 3. The lot number and the expiration date are printed on the labeling.

**STORAGE AND STABILITY**  
 1. Store as packaged in the hermetic bag at the temperature (2-30°C or 38-65°F) and avoid direct sunshine.

**TEST PROCEDURE**  
 1. Sputum sample: Use the cotton swab to pick up 10-50mg sputum samples (equivalent to the size of a match head). Open the cap of sample extraction tube, break the swab tip into the tube, close the disposable sample extraction tube and shake to mix the sample completely. Leave the swab in the extraction tube for one minute.  
 2. Throat swab sampling: saline can be collected in the throat using a larger cotton swab (optional). However, given the greater reliability of sputum sampling, the method described above is recommended.  
 3. Stool sample: Unscrew the sampling bottle, use the sampling rod to pick up 10-50mg of fresh stool samples (equivalent to the size of a match head). Put them into the tube and shake to mix the sample completely.  
 4. Place the test card from the packaging bag, placed on a stable, cut off the protrusion of the collection tube, and add 3 drops of the sample solution into the sample loading hole vertically.  
 5. Read the result in 15 minutes. The result is considered inaccurate and invalid after 30 minutes.

**RESULTS OF INTERPRETATION**

**POSITIVE**  
 Both of C and T lines are appeared in 15 minutes.

**NEGATIVE**  
 C line is appeared while no T line appeared in 15 minutes after the sample is loaded.

**INVALID**  
 As long as the C line does not appear, it indicates that the test result is invalid, and should repeat the sample with another test card.

**NOTE**  
 1. The Novel Coronavirus Antigen Detection Kit (Colloidal Gold) is applicable to stool or sputum samples. Blood, serum, plasma, urine and other samples may cause abnormal results.  
 2. Please make sure that a proper amount of sample is added for testing. Too much or too little sample amount may cause deviations in results.  
 3. For positive judgment, it can be confirmed as soon both T and C line appeared. That may take 3-15 minutes after the sample loading. For the negative judgment, please wait for 15 minutes after sample loading. The result is invalid after 30 minutes after sample loading.  
 4. This product is disposable. DO NOT recycle.  
 5. Dispose of used products, samples, and other consumables as medical wastes under relevant regulations.  
 6. If the test line or control line is out of the test window, do not use the test card. The test result is invalid and repeat the sample with another one.

**INDEX OF SYMBOLS**

<b>IDV</b>	For in vitro diagnostic use		Consult instruction for use
<b>REF</b>	Catalogue number		Temperature limit at 20°C - 30°C
<b>LOT</b>	Batch code		Content sufficient for 25 tests
	Manufacturer		Do not re-use
	Use by		

**FAQ**

**1. What is Coronavirus Antigen test?**  
 Coronavirus Antigen test detects the virus itself by the specific proteins (i.e. Antigens) on its surface, so as to identify the existence of the virus.

**2. What is the difference between Coronavirus Antigen test and Antibody test?**

Antibody test is to detect the antibodies (i.e. antibodies) that people produce after being infected with the virus. Antibody may take several days to produce after infection, and usually positive in the blood for several weeks after the patient recovers. That is, the control panel can also be tested positive for a longer period of time. Therefore, the antibody test is useful in the diagnosis of COVID-19, but with limitations. The Antigen test is to detect the virus, which directly relates to the infection itself.

**3. What is the difference between Antigen test and nucleic acid test?**

Antigen test is based on immune reaction, while nucleic test is based on genetic information. Therefore, Antigen test can be much faster than nucleic acid (PCR) (5 minutes vs 4 hours). Certainly, the cost of PCR is the sensitivity. That is, Antigen test is not sensitive as nucleic acid test. PCR can detect micro amount of SARS-CoV-2 virus due to the amplification process of virus genetic information. This is also why PCR test often takes several hours and may cause aerosol contamination.

**4. Advantages of Antigen test and precautions in application.**

Except for advantages of rapid test, it does not require laboratory processing and medical professionals to operate. Therefore, it is very suitable for large-scale general screening when lacking laboratory conditions and medical professionals. However, if human body contains a very low amount of virus, and with no obvious clinical symptoms, they are likely to be tested negative. The ability of such infected people to spread the virus is also much lower than that of ordinary patients. Therefore, for asymptomatic infected people, it is recommended to adopt nucleic acid for detection. Antigen test helps to quickly identify people with high loads of infection (those who are most likely to infect others). Therefore, Antigen test can help prevent the spread of the pandemic, by identifying those who are most likely to spread the disease and isolating them from others. Thus, shift the focus to identify the most infectious people.

**5. Precautions for Antigen test operations**

Antigen test is to detect the virus itself, so sample selection is very important. We recommend sputum sample, since sputum is the secretion from respiratory tract, whose viruses cluster the most, so sputum is relatively easy and more accurate for detection versus the throat sampling. Because of the uneven distribution of the virus of the respiratory tract mucosa appearance, it is very likely that no virus of very small amount of virus is taken while sampling and lead to false negative result, when sampling via throat swab. Therefore, the accuracy will be lower than of sputum. In addition, studies have shown that other people are infected, the viruses can be detected in the digestive tract, especially in the intestinal stool, therefore if convenient (e.g. in hospital or at home), stool sample is another option for testing. Stool samples mixed with sputum also have a certain probability of detection, although the amount of virus contained is lower than pure sputum. For false negative testing result (PCR positive, Antigen test negative) patient has asymptomatic infection, it is recommended to take nucleic acid as standard, and if patients caught of have other symptoms, it's suggested to re-sampling and test again. Besides, S protein shows good specificity and usually does not cause false positive results.

**6. Others**

The WHO has always been a supporter of Antigen testing. Experts pointed out that compared with conducting an accurate test every two weeks, conducting a relatively insensitive quick test twice a week can more effectively contain the spread of SARS-CoV-2. That is, the focus of the virus detection should be on identifying those who are likely to spread SARS-CoV-2, rather than locating anyone who is infected with SARS-CoV-2. Delhi is the first place in India to start using rapid Antigen testing in June. By mid-July, the number of cases there had been reduced, and the daily death toll had also stabilized, indicating that Antigen testing played a certain role in controlling the spread of the virus. Although after cases declined, India began to lift blockade restrictions and cause the infection to rise again. The Philippine Association of Microbiology and Infectious Diseases has issued interim guidelines for clinicians and medical staff stating that Antigen testing can be used as a substitute for PCR to diagnose coronavirus infection in people with symptoms in the early stage, especially during the acute infection period.

**Clinical Study Report - Abstract**

Home of IDV: Novel Coronavirus Antigen Detection Kit (Colloidal Gold)  
**REF: COVID-19/NG02**  
**Specifications: 25 tests / box**  
**Report Date: August 29th, 2020**

**Summary**

The Novel Coronavirus Antigen Detection Kit (Colloidal Gold) can quickly and qualitatively detect the nucleocapsid protein of novel coronavirus (SARS-CoV-2) in human sputum/stool samples. It can be used as a supplementary test for COVID-19 diagnosis. According to the clinical trial plan, the Novel Coronavirus Antigen Detection Kit (Colloidal Gold) or "test reagent" is to test sputum/stool samples from COVID-19 suspects. Test results are compared with another commercial SARS-CoV-2 nucleic acid detection kit with NMP approval, which is defined as the "gold standard." The sensitivity, specificity, and total agreement rate are used to evaluate the reliability of the test reagent in clinical applications. Method: A collection of clinical samples are examined by the test reagent and the gold standard in parallel, to calculate the clinical sensitivity, clinical specificity and total agreement rate of the test reagent. Standard of criteria for a qualified test reagent: Clinical sensitivity ≥90%, clinical specificity ≥90%, and total agreement rate ≥90%.

**Results**

Compared to the gold standard, the clinical sensitivity of test reagent is 96.0%, the clinical specificity is 97.0%, and the total agreement rate is 96.7%.

**Conclusion**

Compared to the gold standard reagent, the test reagent has reliable performance in diagnosing COVID-19 cases. The sampling method with sputum sample for the respiratory tract mucosa appearance, it is likely that no virus or very small amount of virus is taken while sampling and lead to false negative results, when sampling with throat swab. Therefore the accuracy of throat swab sampling is lower than with sputum sampling.

Distributor:  
 AXELMED SRL  
 Via della Libertazione, 58  
 20098 San Giuliano Milanese (MI)  
 Tel: 0298282694 - info@axelmed.com www.myaxelmed.com  
 Certified Company ISO 13485:2016